ChromoGel™ Candida Agar

The medium use for isolation and identification of Candida bacteria
(Opaque background)

Formula	g/l
Peptone, Yeast extract	24.5
Growth factor	20
Chromogenic mix	5
Chloramphenicol	0.5
Agar	15
Total	65
pH: 6.0 ± 0.2 at 25°C	

Preparation

- # Dissolve 65g in 1000ml of distilled water (purified water)
- # Heat and frequent stirring to completely dissolve the medium
- # Heat with frequent stirring until boiling (**DO NOT AUTOCLAVE**)
- # Cool to 45-50°C, mix well and dispense into appropriate sterile containers

Description

The medium can be used according to usual laboratory techniques with incubation at 35°C-37°C and results reading after 24 and 48 hours. The chromatic characteristics of the colonies after 48 h of incubation at 35°C-37°C are the following:

- ✓ green, blue-green colonies: characteristic of *C.albicans*
- ✓ blue-grey colonies with violet shade and/or violet halo: characteristic of *C.tropicalis*
- ✓ pink, pink-violet, large, rough, spreading, colonies: characteristic of *C.krusei*
- ✓ White or pink or pink-violet colonies: characteristic of other *Candida* species.

Confirm and complete the colonies identification with suitable tests.

Quality Control (37°C-24 H)

Organism	Ref. No.	Result
<i>C. albicans</i>	ATCC 10231	Green - blue colonies
<i>C. tropicalis</i>	NCPF 8841	Blue grey colonies

Technique**Membrane filtration method:**

1. Filter the sample on a 0.45µm membrane and put into surface agar plate.
2. Aerobic Incubate: 18-24h
3. Count *Candida* bacteria

Streaking method:

1. Streaking 0.1-1ml of the sample on surface agar plate.

2. Aerobic Incubate: 18-24h

3. Count *Candida* bacteria

Spreading method:

1. Spreading 0.1-1ml of the sample on surface agar plate.
2. Aerobic Incubate: 18-24h
3. Count *Candida* bacteria

Quality Control acc. to ISO 11133

- Productivity control at 37 ± 2°C / 18-24H

Organism	Inoculated CFU	Result
<i>Candida albicans</i> ATCC 10231	10 ³ - 10 ⁵	Green colonies >30 CFU
<i>Candida tropicalis</i> NCPF 8841	10 ³ - 10 ⁵	Blue grey colonies >30 CFU

- Selective control at 37 ± 2°C / 18-24H

Organism	Inoculated CFU	Result
<i>Enterococcus faecalis</i> ATCC 19433	10 ⁶ - 10 ⁸	Inhibited
<i>E. coli</i> ATCC 25922	10 ⁶ - 10 ⁸	Inhibited


Storage

Store the medium tightly closed in a dry place at room temperature and use before the expiry date on the label.

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