ChromoGel™ E.coli Coliforms Agar – ECC

ChromoGel™ E.coli Coliforms Agar is a selective medium used for the isolation and identification of E.coli and Coliforms from clinical specimens, foodstuff and environmental samples.

| Formula | g/l |
|------------------------------|-----------|
| Peptone & Yeast extract | 18 |
| Growth factor | 22 |
| Chromogenic & Selective Mix | 3 |
| Agar | 15 |
| Total: | 58 |
| pH: 7.2 ± 0.2 at 25°C | |

Preparation

- # Dissolve 58g in 1000ml of distilled water (purified water)
- # Heat and frequent stirring to completely dissolve the medium
- # Autoclave at 121°C for 15 minutes
- # Mix well and dispense into appropriate sterile containers

Description

ChromoGel™ E coli Coliforms Agar is a selective medium used for the isolation and identification of E.coli and Coliforms from clinical specimens, foodstuff and environmental samples.

- The medium can use as secondary selective agar plate on ISO method, petri-film method.
- The samples can streak directly on surface agar plates without enrichment step.
- Aerobic Incubate:
 - At 30 ± 2°C for 18-24 hours for Total Coliforms
 - At 37 ± 2°C for 18-24 hours for E.coli and Coliforms
 - At 44 ± 2°C for 18-24 hours for faecal Coliforms
- Read the result within 18-24h

| Bacteria strain | Chromogenic mix | Colonies Color |
|-----------------|-----------------|-------------------------------|
| E.coli | + | Purple |
| Coliforms | + | Blue to green |
| Other bacteria | - | Colourless White or inhibited |

Technique**Membrane filtration method:**

1. Filter the sample on a 0.45µm membrane and put into surface agar plate.
2. Aerobic Incubate: 18-24h
3. Count as E. coli and coliforms confirmed by indole test.

Streaking method:

1. Streaking 1-10µl of the sample on surface agar plate.

2. Aerobic Incubate: 18-24h
3. Count as E. coli and coliforms confirmed by indole test.

Spreading method:

1. Spreading 0.1-1ml of the sample on surface agar plate.
2. Aerobic Incubate: 18-24h
3. Count as E. coli and coliforms confirmed by indole test.

Quality Control acc. to ISO 11133

- Productivity control at 37 ± 2°C / 18-24H

| Organism | Inoculated CFU | Result |
|--------------------------------------|-----------------------------------|--------------------------|
| E.coli ATCC 25922 | 10 ³ - 10 ⁵ | Purple >30 CFU |
| E.aerogenes ATCC 13048 | 10 ³ - 10 ⁵ | Blue to green >30 CFU |
| Klebsiella pneumoniae NCTC 9633 | 10 ³ - 10 ⁵ | Blue to green >30 CFU |
| Salmonella typhimurium NCTC 12023 | 10 ³ - 10 ⁵ | White >30 CFU |

- Selective control at 37 ± 2°C / 18-24H

| Organism | Inoculated CFU | Result |
|----------------------------|-----------------------------------|-----------|
| E. faecalis ATCC 19433 | 10 ⁶ - 10 ⁸ | Inhibited |
| Staph aureus NCTC 12981 | 10 ⁶ - 10 ⁸ | Inhibited |

Storage

Dehydrate culture media powder keep at 2°C to 30 °C

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