

**AgarCult™
Egg Yolk Emulsion**

Formula	100ml
Egg Yolk	50ml
Sodium chloride solution	50ml

Application medium

AgarCult™ Bacillus Cereus Agar – MYP

AgarCult™ Bacillus Cereus Agar - PEMBA

AgarCult™ TSC Agar

Preparation

For MYP Medium: Fill 1 vial supplement into 450ml of sterile MYP medium and 50ml of Egg Yolk Emulsion at 45-50°C. Mix and dispense into sterile containers.

For PEMBA Medium: Fill 1 vial supplement into 475ml of sterile PEMBA medium and 25ml of Egg Yolk Emulsion at 45-50°C. Mix and dispense into sterile containers.

For TSC Medium: Fill 1 vial supplement into 450ml of sterile TSC medium and 50ml of Egg Yolk Emulsion at 45-50°C. Mix and dispense into sterile containers.

Storage

Store the Supplement in dry place at 2-8°C and use before the expiry date on the label.

Ordering Information

Code	Products	Package
EYE50	AgarCult™ Egg Yolk Emulsion, 50ml	Bottle
EYE100	AgarCult™ Egg Yolk Emulsion, 100ml	Bottle
EYE10100	AgarCult™ Egg Yolk Emulsion, 10x100ml	Bottle