

AgarCult™
Egg Yolk Tellurite

Formula	100ml
Egg Yolk	20ml
Potassium tellurite	0.25g
Sodium chloride	0.75g
Distilled water	80ml

Application medium

AgarCult™ Baird Parker Agar Base

Preparation

Fill 25ml of Egg Yolk Tellurite into 475ml Bair Parker medium at 45-50°C.
Mix and dispense into sterile containers.

Storage

Store the Supplement in dry place at 2-8°C and use before the expiry date on the label.

Ordering Information

Code	Products	Package
EYT50	AgarCult™ Egg Yolk Tellurite, 50ml	Bottle
EYT100	AgarCult™ Egg Yolk Tellurite, 100ml	Bottle
EYT10100	AgarCult™ Egg Yolk Tellurite, 10x100ml	Bottle