

**AgarCult™**  
**Egg Yolk Emulsion with Polymycin B**

Formula	100ml
Egg Yolk	50ml
Polymycin B	100.000UI
Sodium chloride solution	50ml

**Application medium**

AgarCult™ Bacillus Cereus Agar – MYP

AgarCult™ Bacillus Cereus Agar - PEMBA

**Preparation**

# For MYP Medium: Mix 450ml of sterile MYP medium and 50ml of Egg Yolk Emulsion with Polymycin B at 45-50°C. Mix and dispense into sterile containers.

# For PEMBA Medium: Mix 475ml of sterile PEMBA medium and 25ml of Egg Yolk Emulsion with Polymycin B at 45-50°C. Mix and dispense into sterile containers.

**Storage**

Store the Supplement in dry place at 2-8°C and use before the expiry date on the label.

**Ordering Information**

Code	Products	Package
EYEP50	AgarCult™ Egg Yolk Emulsion w Polymycin B, 50ml	Bottle
EYEP100	AgarCult™ Egg Yolk Emulsion w Polymycin B, 100ml	Bottle
EYEP10100	AgarCult™ Egg Yolk Emulsion w Polymycin B, 10x100ml	Bottle